

Pulmicort Flexhaler

Pulmicort Flexhaler®

90 mcg (green on label)

180mcg (purple label))



Medicine	Need a	My dose
Controller (dry powder inhaler)	No	Check Asthma Action Plan



Figure 2 - Load



Figure 3 - Twist



Figure 4 - Click




Figure 5 - Inhale

1. Holding Flexhaler upright, remove cap
2. Twist the brown grip until you hear a click and counter has decreased by 1, then click back
3. Place mouth over mouthpiece with a tight seal, holding device parallel to the ground
4. Inhale quickly and deeply
5. Hold breath as long as you can, around 10 seconds
6. Remove Flexhaler from mouth and exhale slowly away from device.
7. Place cap back on and twist until you hear a click and arrow is lined up with the counter
8. Repeat steps 1-7 for each dose



Use the camera on your phone to scan this QR code for the Pulmicort Flexhaler Website:

There is information about:

- ◆ How to use the Pulmicort Flexhaler